



- ✓ Neutral colouring
- ✓ Does not inhibit adhesion of other materials
- ✓ Shows clearly on X-rays



AHseal

Powder and resin root canal filling material which has excellent sealing capacity. It offers close adaptations of the root canal walls and has very low contraction during setting.

Setting time 8 - 20 hours.

Working time over 4 hours.

Indications

Permanent obturation of root canals of teeth of the secondary dentition with or without the aid of obturation points.

Features

Long working time during preparation
Radiopaque

Advantages

Neutral colouring
Does not inhibit adhesion of other materials
Aids material manipulation outside mouth
Shows clearly on X-rays

Available in	Code
10gm Resin/8gm Powder	AH0780

Your global partner

Delivering quality for more than 30 years

Advanced Healthcare Ltd.

Units 2-4 Leavers Estate, Chiddingstone Causeway, Tonbridge, Kent TN11 8JU UK

Tel: +44(0)1892 870500 Fax: +44(0)1892 870482 Email: sales@ahl.uk.com

www.ahl.uk.com

Authorised distributor

© 2023 Advanced Healthcare Ltd. E&OE.