





✓ Good early translucency

High strength

- Packable consistency
- Good wear and erosion resistance
- ✓ Vita® shades
- Easy mixing

Powder/Liquid 15gm/7ml plus mixing pad and scoop









AHfil Anterior Restorative Material

AHfil glass ionomer filling materials are manufactured using translucent aluminasilicate glass that offers good aesthetics without compromising on strength and wear resistance.

AHfil anterior restorative material delivers good translucency and is biocompatible with minimal pulp sensitivity. It has a CTE close to that of natural dentition and good bond strength. Optimum working and setting time. Easy mixing and handling.

Indications

Class V cavities
Class III cavities
Class I and II cavities in
deciduous teeth
Pits and fissures

Features

Good early translucency High strength Packable consistency Good wear and erosion resistance 5 VITA® shades Easy mixing

Advantages

Good aesthetics Long lasting restorations Ease of handling Durable restorations Good tooth matching

Code
AH0001
AH0002
AH0003
AH0004

Property	ISO9917-1:2007 Classes B & C	Typical Properties
Working Time	_	1 min 57 sec
Net Setting Time	1 min 30 secs – 6 mins	2 mins 32 secs
Compressive Strength	100MPa Minimum	200+ MPa

Your global partner

Delivering quality for more than 30 years

Advanced Healthcare Ltd.

Units 2-4 Leavers Estate, Chiddingstone Causeway, Tonbridge, Kent TN11 8JU UK Tel: +44(0)1892 870500 Fax: +44(0)1892 870482 Email: sales@ahl.uk.com