



5ml syringe with 10 Luer lock application tips

caviEtch Phosphoric Acid Etchant

caviEtch is a 40% aqueous phosphoric acid etchant designed to deliver a high quality "etch" to enamel surfaces for bonding to composite and resin restorations.

The etch gel produces micromechanical interlock for the keying of resin based cements and restoratives. The gel is thixotropic meaning it can be accurately placed on the area to be treated without risk of spreading. The bright green colour clearly shows the clinician where the etchant is placed and also useful when it comes to washing off the product.

Optimum viscosity for easy extrusion with thixotropic non-drip effect.
Indicated for condition of enamel and dentin prior to adhesive materials.
Bright green colour for precise placement.

Indications

Condition of enamel (and dentin where the total etch technique is desired) prior to adhesive materials placement

Features and benefits

Treatment of enamel margins of composite restorations
Treatment of enamel surfaces of composite restorations
Treatment of occlusal surface prior to application of sealant resin
Treatment of enamel surfaces of prior to the attachment of orthodontic appliances by means of resin adhesive

Available in	Code
5ml syringe	AH1300

Your global partner

Delivering quality for more than 30 years

Advanced Healthcare Ltd.

Units 2-4 Leavers Estate, Chiddingstone Causeway, Tonbridge, Kent TN11 8JU UK

Tel: +44(0)1892 870500 Fax: +44(0)1892 870482 Email: sales@ahl.uk.com

www.ahl.uk.com

Authorised distributor

© 2023 Advanced Healthcare Ltd. E&OE.